

## 兔抗 NSF 多克隆抗体

中文名称：兔抗 NSF 多克隆抗体

英文名称：Anti-NSF rabbit polyclonal antibody

别名：SKD2

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：NSF

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

<b>Background:</b>	Required for vesicle-mediated transport. Catalyzes the fusion of transport vesicles within the Golgi cisternae. Is also required for transport from the endoplasmic reticulum to the Golgi stack. Seems to function as a fusion protein required for the delivery of cargo proteins to all compartments of the Golgi stack independent of vesicle origin. Interaction with AMPAR subunit GRIA2 leads to influence GRIA2 membrane cycling By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NSF
<b>Immunogen:</b>	Synthesized peptide derived from internal of human NSF.
<b>Full name:</b>	N-ethylmaleimide-sensitive factor
<b>Synonyms :</b>	SKD2

<b>SwissProt:</b>	P46459
<b>WB Predicted band size:</b>	83 kDa
<b>WB Positive control:</b>	COLO cells lysate
<b>WB Recommended dilution:</b>	500-3000

