

## 兔抗 IL18R1 多克隆抗体

- 中文名称：兔抗 IL18R1 多克隆抗体
- 英文名称：Anti-IL18R1 rabbit polyclonal antibody
- 别名：CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp
- 抗原：IL18R1
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

**Background:**

The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor specifically binds interleukin 18 (IL18), and is essential for IL18 mediated signal transduction. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. This gene along with four other members of the interleukin 1 receptor family, including IL1R2, IL1R1, ILRL2 (IL-1Rrp2), and IL1RL1 (T1/ST2), form a gene cluster on chromosome 2q. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IL18R1
<b>Immunogen:</b>	Fusion protein of human IL18R1
<b>Full name:</b>	interleukin 18 receptor 1
<b>Synonyms :</b>	CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp
<b>SwissProt:</b>	Q13478
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>IHC Recommend dilution:</b>	50-200

