

## 兔抗 RPL10L 多克隆抗体

中文名称：兔抗 RPL10L 多克隆抗体

英文名称：Anti-RPL10L rabbit polyclonal antibody

别名：ribosomal protein L10 like; RPL10\_5\_1358

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RPL10L

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein sharing sequence similarity with ribosomal protein L10. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. This gene is intronless.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RPL10L
<b>Immunogen:</b>	Fusion protein of human RPL10L
<b>Full name:</b>	ribosomal protein L10 like
<b>Synonyms:</b>	RPL10_5_1358
<b>SwissProt:</b>	Q96L21

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-200

