

## 兔抗 NSMCE3 多克隆抗体

中文名称：兔抗 NSMCE3 多克隆抗体

英文名称：Anti-NSMCE3 rabbit polyclonal antibody

别名：NSE3 homolog, SMC5-SMC6 complex component; HCA4; L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NSMCE3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is part of the SMC5-6 chromatin reorganizing complex and is a member of the MAGE superfamily. This is an intronless gene. [provided by RefSeq, May 2011]
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NSMCE3
<b>Immunogen:</b>	Synthetic peptide of human NSMCE3
<b>Full name:</b>	NSE3 homolog, SMC5-SMC6 complex component
<b>Synonyms:</b>	HCA4; LICS; NSE3; NDNL2; MAGEG1; MAGEL3
<b>SwissProt:</b>	Q96MG7
<b>IHC positive control:</b>	Human breast cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200

<b>WB Predicted band size:</b>	34 kDa
<b>WB Positive control:</b>	Hela cells and human prostate tissue, A172 cells and mouse liver tissue
<b>WB Recommended dilution:</b>	1000-5000

