

NADK2 抗原（重组蛋白）

中文名称：NADK2 抗原（重组蛋白）

英文名称：NADK2 Antigen (Recombinant Protein)

别名：NAD kinase 2, mitochondrial; DECRD; MNADK; NADKD1; C5orf33

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 164-442 amino acids of human NADK2

技术规格

Full name:	NAD kinase 2, mitochondrial
Synonyms:	DECRD; MNADK; NADKD1; C5orf33
Swissprot:	Q4G0N4
Gene Accession:	BC062567
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a mitochondrial kinase that catalyzes the phosphorylation of NAD to yield NADP. Mutations in this gene result in 2, 4-dienoyl-CoA reductase deficiency. Alternative splicing results in multiple transcript variants.

mlbio 梅联生物
Good elisakit producers