

兔抗 NT5E 多克隆抗体

中文名称：兔抗 NT5E 多克隆抗体

英文名称： Anti-NT5E rabbit polyclonal antibody

别名： NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA

相关类别： 一抗

抗原： NT5E

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	NT5E
Immunogen:	Fusion protein of human NT5E
Full name:	5'-nucleotidase, ecto (CD73)
Synonyms :	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA

SwissProt:	P21589
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	63 kDa
WB Positive control:	A549 cells
WB Recommended dilution:	500-2000



