

NAT10 抗原(重组蛋白)

中文名称: NAT10 抗原(重组蛋白)

英文名称: NAT10 Antigen (Recombinant Protein)

别 名: N-acetyltransferase 10 (GCN5-related); ALP; NET43

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 558-753 amino acids of human NAT10

技术规格

Full name:	N-acetyltransferase 10 (GCN5-related)
Synonyms:	ALP; NET43
Swissprot:	Q9H0A0
Gene Accession:	BC035558
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	N-acetyltransferase 10 is an enzyme that in humans is encoded by the NAT10 gene. Has protein acetyltransferase activity in vitro. Can acetylate both histones and microtubules. Histone acetylation may regulate transcription and mitotic chromosome de-condensation. Ac tivates telomerase activity by stimulating the transcription of TERT, and may also regulate telomerase function by affecting the balance of telomerase subunit assembly, disassembly, and localization.