

## 兔抗 IL20 多克隆抗体

中文名称：兔抗 IL20 多克隆抗体

英文名称：Anti-IL20 rabbit polyclonal antibody

别名：interleukin 20; IL-20; IL10D; ZCYTO10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IL20

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	IL20
<b>Immunogen:</b>	Synthetic peptide of human IL20
<b>Full name:</b>	interleukin 20
<b>Synonyms:</b>	IL-20; IL10D; ZCYTO10
<b>SwissProt:</b>	Q9NYY1
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	25-100

