

NCLN 抗原（重组蛋白）

中文名称：NCLN 抗原（重组蛋白）

英文名称：NCLN Antigen (Recombinant Protein)

别名：nicalin; NET59

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 364-563 amino acids of human NCLN

技术规格

Full name:	nicalin
Synonyms:	NET59
Swissprot:	Q969V3
Gene Accession:	BC013283
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May antagonize Nodal signaling and subsequent organization of axial structures during mesodermal patterning.