

兔抗 IL23R 多克隆抗体

中文名称:	兔抗 IL23R 多克隆抗体
英文名称:	Anti-IL23R rabbit polyclonal antibody
储存:	冷冻 (-20℃)
抗原:	IL23R
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a subunit of the receptor for IL23A/IL23. This protein pairs with the receptor molecule IL12RB1/IL12Rbeta1, and both are required for IL23A signaling. This protein associates constitutively with Janus kinase 2 (JAK2), and also binds to transcription activator STAT3 in a ligand-dependent manner.
Applications:	ELISA, IHC
Name of antibody:	IL23R
Immunogen:	Fusion protein of human IL23R
Full name:	Interleukin-23 receptor
SwissProt:	Q5VWK5
ELISA Recommended dilution:	500-5000

IHC positive control:	Human lung cancer
IHC Recommend dilution:	10-50

