

NCSTN 抗原（重组蛋白）

中文名称： NCSTN 抗原（重组蛋白）

英文名称： NCSTN Antigen (Recombinant Protein)

别名： nicastrin; ATAG1874

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 34-330 amino acids of human NCSTN

技术规格

Full name:	nicastrin
Synonyms:	ATAG1874
Swissprot:	Q92542
Gene Accession:	BC047621
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a type I transmembrane glycoprotein that is an integral component of the multimeric gamma-secretase complex. The encoded protein cleaves integral membrane proteins, including Notch receptors and beta-amyloid precursor protein, and may be a stabilizing cofactor required for gamma-secretase complex assembly.