

兔抗 NTM 多克隆抗体

中文名称: 兔抗 NTM 多克隆抗体

英文名称: Anti-NTM rabbit polyclonal antibody

别名: HNT; NTRI; IGLON2

储存: 冷冻 (-20℃) 避光

抗原: NTM

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a member of the IgLON (LAMP, OBCAM , Ntm) family of immunoglobulin (Ig) domain-containing glycosylphosphatidylinositol (GPI)-anchored cell adhesion molecules. The encoded protein may promote neurite outgrowth and adhesion via a homophilic mechanism. This gene is closely linked to a related family member, opioid binding protein/cell adhesion molecule-like (OPCML), on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2009]
Applications:	WB
Name of antibody:	NTM
Immunogen:	Synthesized peptide derived from internal of human NT.

Full name:	neurotrimin
Synonyms :	HNT; NTRI; IGLON2
SwissProt:	Q9P121
WB Predicted band size:	38 kDa
WB Positive control:	Jurkat cells, Huvec cells and HepG2 cells lysates
WB Recommended dilution:	500-3000

