

兔抗 RPL26L1 多克隆抗体

中文名称：兔抗 RPL26L1 多克隆抗体

英文名称： Anti-RPL26L1 rabbit polyclonal antibody

别名： RPL26P1

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： RPL26L1

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. Transcript variants utilizing alternative polyA signals exist.
Applications:	ELISA, IHC
Name of antibody:	RPL26L1
Immunogen:	Fusion protein of human RPL26L1
Full name:	ribosomal protein L26-like 1
Synonyms :	RPL26P1

SwissProt:	Q9UNX3
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and human lung cancer
IHC Recommend dilution:	50-200

