

## 兔抗 RPL27A 多克隆抗体

中文名称：兔抗 RPL27A 多克隆抗体

英文名称： Anti-RPL27A rabbit polyclonal antibody

别名： L27A

抗原： RPL27A

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

**Background:**

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L15P family of ribosomal proteins. It is located in the cytoplasm. Variable expression of this gene in colorectal cancers compared to adjacent normal tissues has been observed, although no correlation between the level of expression and the severity of the disease has been found. As is typical for

	genes encoding ribosomal proteins, multiple processed ps eudogenes derived from this gene are dispersed through the genome.
<b>Applications:</b>	WB, IF
<b>Name of antibody:</b>	RPL27A
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human R PL27A.
<b>Full name:</b>	ibosomal protein L27a
<b>Synonyms :</b>	L27A
<b>SwissProt:</b>	P46776
<b>WB Predicted band size:</b>	17 kDa
<b>WB Positive control:</b>	COLO cells lysate
<b>WB Recommended dilution:</b>	500-3000
<b>IF positive control:</b>	A549 cells
<b>IF Recommend dilution:</b>	100-500



