

兔抗 IL26 多克隆抗体

- 中文名称： 兔抗 IL26 多克隆抗体
- 英文名称： Anti-IL26 rabbit polyclonal antibody
- 别名： interleukin 26; AK155; IL-26
- 相关类别： 一抗
- 抗原： IL26
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 反应种属： Human
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

技术规格

Background:	This gene was identified by its overexpression specifically in herpesvirus samimiri-transformed T cells. The encoded protein is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus samimiri.
Applications:	ELISA, IHC
Name of antibody:	IL26
Immunogen:	Synthetic peptide of human IL26
Full name:	interleukin 26
Synonyms:	AK155; IL-26

SwissProt:	Q9NPH9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150

