

兔抗 UNKL 多克隆抗体

中文名称：兔抗 UNKL 多克隆抗体

英文名称：Anti-UNKL rabbit polyclonal antibody

别名：unkempt family like zinc finger; ZC3H5L; C16orf28; ZC3HDC5L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UNKL

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a RING finger protein that may function in Rac signaling. It can bind to Brg/Brm-associated factor 60b and can promote its ubiquitination. Alternative splicing of this gene results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	UNKL
Immunogen:	Fusion protein of human UNKL
Full name:	unkempt family like zinc finger
Synonyms:	ZC3H5L; C16orf28; ZC3HDC5L
SwissProt:	Q9H9P5

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human esophagus cancer
IHC Recommend dilution:	50-300



