

兔抗 RPL28 多克隆抗体

中文名称：兔抗 RPL28 多克隆抗体

英文名称：Anti-RPL28 rabbit polyclonal antibody

别名：ribosomal protein L28; L28

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RPL28

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L28E family of ribosomal proteins. It is located in the cytoplasm. Variable expression of this gene in colorectal cancers compared to adjacent normal tissues has been observed, although no correlation between the level of expression and the severity of the disease has been found. As is typical fo

	r genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
Applications:	ELISA, IHC
Name of antibody:	RPL28
Immunogen:	Fusion protein of human RPL28
Full name:	ribosomal protein L28
Synonyms:	L28
SwissProt:	P46779
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	40-200

