

兔抗 RPL30 多克隆抗体

中文名称：兔抗 RPL30 多克隆抗体

英文名称： Anti-RPL30 rabbit polyclonal antibody

别名： L30

抗原： RPL30

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L30E family of ribosomal proteins. It is located in the cytoplasm. This gene is co-transcribed with the U72 small nucleolar RNA gene, which is located in its fourth intron. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene d

	ispersed through the genome.
Applications:	WB, IHC
Name of antibody:	RPL30
Immunogen:	Synthesized peptide derived from internal of human RPL30.
Full name:	ribosomal protein L30
Synonyms :	L30
SwissProt:	P62888
IHC positive control:	Human prostate carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	13 kDa
WB Positive control:	A549 cells, COS-7 cells and RAW264.7 cells lysates
WB Recommended dilution:	500-3000



