

## 兔抗 IL27 多克隆抗体

- 中文名称：兔抗 IL27 多克隆抗体
- 英文名称：Anti-IL27 rabbit polyclonal antibody
- 别名：p28; IL30; IL-27; IL27A; IL-27A; IL27p28
- 抗原：IL27
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is one of the subunits of a heterodimeric cytokine complex. This protein is related to interleukin 12A (IL12A). It interacts with Epstein-Barr virus induced gene 3 (EBI3), a protein similar to interleukin 12B (IL12B), and forms a complex that has been shown to drive rapid expansion of naive but not memory CD4(+) T cells. The complex is also found to synergize strongly with interleukin 12 to trigger interferon gamma (IFNG) production of naive CD4(+) T cells. The biological effects of this cytokine are mediated by the class I cytokine receptor (WSX1/TCRR)
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	IL27
<b>Immunogen:</b>	Fusion protein of human IL27
<b>Full name:</b>	interleukin 27
<b>Synonyms :</b>	p28; IL30; IL-27; IL27A; IL-27A; IL27p28
<b>SwissProt:</b>	Q8NEV9
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human lymphoma and human liver cancer
<b>IHC Recommend dilution:</b>	25-100

