

兔抗 RPL31 多克隆抗体

中文名称：兔抗 RPL31 多克隆抗体

英文名称： Anti-RPL31 rabbit polyclonal antibody

别名： L31

抗原： RPL31

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L31E family of ribosomal proteins. It is located in the cytoplasm. Higher levels of expression of this gene in familial adenomatous polyposis compared to matched normal tissues have been observed. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudo

	genes of this gene dispersed through the genome. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	IHC
Name of antibody:	RPL31
Immunogen:	Synthesized peptide derived from internal of human RPL31.
Full name:	ribosomal protein L31
Synonyms :	L31
SwissProt:	P62899
IHC positive control:	Human prostate tissue
IHC Recommend dilution:	50-100

