

## 兔抗 RPL35A 多克隆抗体

中文名称：兔抗 RPL35A 多克隆抗体

英文名称： Anti-RPL35A rabbit polyclonal antibody

别名： DBA5; L35A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPL35A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L35AE family of ribosomal proteins. It is located in the cytoplasm. The rat protein has been shown to bind to both initiator and elongator tRNAs, and thus, it is located at the P site, or P and A sites, of the ribosome.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RPL35A
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	ribosomal protein L35a
<b>Synonyms :</b>	DBA5; L35A
<b>SwissProt:</b>	P18077
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

