

兔抗 RPL35 多克隆抗体

中文名称：兔抗 RPL35 多克隆抗体

英文名称： Anti-RPL35 rabbit polyclonal antibody

别名： ribosomal protein L35; L35

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPL35

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L29P family of ribosomal proteins. It is located in the cytoplasm. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed throughout the genome.
Applications:	ELISA, WB, IHC

Name of antibody:	RPL35
Immunogen:	Fusion protein of human RPL35
Full name:	ribosomal protein L35
Synonyms:	L35
SwissProt:	P42766
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	15 kDa
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	500-2000



