

## 兔抗 UQCC2 多克隆抗体

中文名称：兔抗 UQCC2 多克隆抗体

英文名称：Anti-UQCC2 rabbit polyclonal antibody

别名：ubiquinol-cytochrome c reductase complex assembly factor 2; M19; Cbp6; MNF1; C6orf125; C6orf126; bA6B20.2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UQCC2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclona

### 技术规格

<b>Background:</b>	This gene encodes a nucleoid protein localized to the mitochondria inner membrane. The encoded protein affects regulation of insulin secretion, mitochondrial ATP production, and myogenesis through modulation of mitochondrial respiratory chain activity.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UQCC2
<b>Immunogen:</b>	Fusion protein of human UQCC2
<b>Full name:</b>	ubiquinol-cytochrome c reductase complex assembly factor 2

<b>Synonyms:</b>	M19; Cbp6; MNF1; C6orf125; C6orf126; bA6B20.2
<b>SwissProt:</b>	Q9BRT2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	40-200

