

兔抗 RPL39L 多克隆抗体

中文名称: 兔抗 RPL39L 多克隆抗体

英文名称: Anti-RPL39L rabbit polyclonal antibody

别名: RPL39L1

储存: 冷冻 (-20℃) 避光

抗原: RPL39L

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a protein sharing high sequence similarity with ribosomal protein L39. Although the name of this gene has been referred to as 'ribosomal protein L39' in the public databases, its official name is 'ribosomal protein L39-like'. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome.
Applications:	IHC
Name of antibody:	RPL39L
Immunogen:	Synthesized peptide derived from internal of human RL39L.
Full name:	ribosomal protein L39-like

Synonyms :	RPL39L1
SwissProt:	Q96EH5
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

