

兔抗 IL36B 多克隆抗体

中文名称: 兔抗 IL36B 多克隆抗体

英文名称: Anti-IL36B rabbit polyclonal antibody

别名: interleukin 36 beta; FIL1; FIL1H; IL1F8; IL1H2; IL-1F8; IL-1H2; IL1-ETA; FIL1-(ETA); FILI-(ETA)

储存: 冷冻 (-20℃)

抗原: IL36B

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the interleukin 1 cytokine family. Protein structure modeling indicated that this cytokine may contain a 12-stranded beta-trefoil structure that is conserved between IL1A (IL-A alpha) and IL1B (IL-1 beta). This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. Two alternatively spliced transcript variants encoding distinct isoforms have been reported.
Applications:	ELISA, IHC
Name of antibody:	IL36B

Immunogen:	Synthetic peptide of human IL36B
Full name:	interleukin 36 beta
Synonyms:	FIL1; FIL1H; IL1F8; IL1H2; IL-1F8; IL-1H2; IL1-ETA; FIL1-(ETA); FILI-(ETA)
SwissProt:	Q9NZH7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

