

## 兔抗 UQCR10 多克隆抗体

中文名称：兔抗 UQCR10 多克隆抗体

英文名称：Anti-UQCR10 rabbit polyclonal antibody

别名：ubiquinol-cytochrome c reductase, complex III subunit X QCR9; UCRC; HSPC051; HSPC119; HSPC151; UCCR7.2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UQCR10

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclona

### 技术规格

<b>Background:</b>	UCRC is a subunit of mitochondrial complex III (ubiquinol-cytochrome c reductase; EC 1.10.2.2), which forms the middle segment of the respiratory chain of the inner mitochondrial membrane (Schagger et al., 1995 [PubMed 8592474]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UQCR10
<b>Immunogen:</b>	Synthetic peptide of human UQCR10
<b>Full name:</b>	ubiquinol-cytochrome c reductase, complex III subunit X
<b>Synonyms:</b>	QCR9; UCRC; HSPC051; HSPC119; HSPC151; UCCR7.2

<b>SwissProt:</b>	Q9UDW1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	20-100



