

## 兔抗 RPLP2 多克隆抗体

中文名称：兔抗 RPLP2 多克隆抗体

英文名称： Anti-RPLP2 rabbit polyclonal antibody

别名： P2; LP2; RPP2; D11S2243E

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPLP2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal phosphoprotein that is a component of the 60S subunit. The protein, which is a functional equivalent of the E. coli L7/L12 ribosomal protein, belongs to the L12P family of ribosomal proteins. It plays an important role in the elongation step of protein synthesis. Unlike most ribosomal proteins,

	ns, which are basic, the encoded protein is acidic.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RPLP2
<b>Immunogen:</b>	Fusion protein of human RPLP2
<b>Full name:</b>	ribosomal protein, large, P2
<b>Synonyms :</b>	P2; LP2; RPP2; D11S2243E
<b>SwissProt:</b>	P05387
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	12 kDa
<b>WB Positive control:</b>	Raji, K562, A549, hepg2, PC3 and hela cells
<b>WB Recommended dilution:</b>	1000-5000



