

兔抗 NUB1 多克隆抗体

中文名称：兔抗 NUB1 多克隆抗体

英文名称： Anti-NUB1 rabbit polyclonal antibody

别名： BS4; NUB1L; NYREN18

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： NUB1

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that functions as a negative regulator of NEDD8, a ubiquitin-like protein that conjugates with cullin family members in order to regulate vital biological events. The protein encoded by this gene regulates the NEDD8 conjugation system post-transcriptionally by recruiting NEDD8 and its conjugates to the proteasome for degradation.
Applications:	WB, IHC
Name of antibody:	NUB1
Immunogen:	Synthesized peptide derived from C-terminal of human NYREN18.
Full name:	negative regulator of ubiquitin-like proteins 1
Synonyms :	BS4; NUB1L; NYREN18

SwissProt:	Q9Y5A7
IHC positive control:	Human tonsil tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	71 kDa
WB Positive control:	COS7 cells lysate
WB Recommended dilution:	500-3000



