

兔抗 RPN1 多克隆抗体

中文名称：兔抗 RPN1 多克隆抗体

英文名称： Anti-RPN1 rabbit polyclonal antibody

别名： OST1; RBPH1

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： RPN1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a type I integral membrane protein found only in the rough endoplasmic reticulum. The encoded protein is part of an N-oligosaccharyl transferase complex that links high mannose oligosaccharides to asparagine residues found in the Asn-X-Ser/Thr consensus motif of nascent polypeptide chains. This protein forms part of the regulatory subunit of the 26S proteasome and may mediate binding of ubiquitin-like domains to this proteasome.
Applications:	ELISA, WB, IHC
Name of antibody:	RPN1

Immunogen:	Fusion protein of human RPN1
Full name:	ribophorin I
Synonyms :	OST1; RBPH1
SwissProt:	P04843
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	69 kDa
WB Positive control:	Hela and A549 cells, human liver cancer tissue
WB Recommended dilution:	500-2000



