

兔抗 NUBP2 多克隆抗体

中文名称：兔抗 NUBP2 多克隆抗体

英文名称： Anti-NUBP2 rabbit polyclonal antibody

别名： CFD1; NBP 2; NUBP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NUBP2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an adenosine triphosphate (ATP) and metal-binding protein that is required for the assembly of cytosolic iron-sulfur proteins. The encoded protein functions in a heterotetramer with nucleotide-binding protein 1 (NUBP1). Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	NUBP2
Immunogen:	Full length fusion protein
Full name:	nucleotide binding protein 2
Synonyms:	CFD1; NBP 2; NUBP1
SwissProt:	Q9Y5Y2
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human ovarian cancer;Human gastric cancer
IHC Recommend dilution:	25-100



