

## 兔抗 ILK(Ab-246) 多克隆抗体

- 中文名称：兔抗 ILK(Ab-246) 多克隆抗体
- 英文名称：Anti-ILK(Ab-246) rabbit polyclonal antibody
- 别名：P59; ILK-1; ILK-2; p59ILK; HEL-S-28
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：ILK(Ab-246)
- 反应种属：Human Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein with a kinase-like domain and four ankyrin-like repeats. The encoded protein associates at the cell membrane with the cytoplasmic domain of beta integrins, where it regulates integrin-mediated signal transduction. Activity of this protein is important in the epithelial to mesenchymal transition, and over-expression of this gene is implicated in tumor growth and metastasis. Alternative splicing results in multiple transcript variants
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	ILK(Ab-246)
<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human ILK aro

	und the phosphorylation site of serine 246 (I-F-S(p)-H-P).
<b>Full name:</b>	integrin-linked kinase
<b>Synonyms :</b>	P59; ILK-1; ILK-2; p59ILK; HEL-S-28
<b>SwissProt:</b>	Q13418
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	51 kDa
<b>WB Positive control:</b>	Rat liver tissue lysate
<b>WB Recommended dilution:</b>	500-3000



