

## 兔抗 URB1 多克隆抗体

中文名称：兔抗 URB1 多克隆抗体

英文名称：Anti-URB1 rabbit polyclonal antibody

级别：NPA1; C21orf108

储存：冷冻 (-20℃)

抗原：URB1

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	URB1, also named as C21orf108, KIAA0539 and NPA1, is a 2,271 amino acid protein that is homologous to the <i>S. cerevisiae</i> protein URB1 which is a nucleolar protein involved in the normal accumulation of 25S and 5.8S rRNAs and may also play a role in the biosynthesis of the 60S ribosomal unit. Also localized to the nucleolus, human URB1 is thought to have the same functions as <i>S. cerevisiae</i> URB1. The antibody is specific to URB1.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	URB1
<b>Immunogen:</b>	Synthetic peptide of human URB1
<b>Full name:</b>	URB1 ribosome biogenesis homolog
<b>Synonyms:</b>	NPA1; C21orf108

<b>SwissProt:</b>	O60287
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human tonsil and Human ovarian cancer
<b>IHC Recommend dilution:</b>	20-100



