

兔抗 URGCP 多克隆抗体

中文名称: 兔抗 URGCP 多克隆抗体

英文名称: Anti-URGCP rabbit polyclonal antibody

别名: URG4

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: URGCP

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	URG4 is upregulated in the presence of hepatitis B virus (HBV)-encoded X antigen (HBxAg) and may contribute to the development of hepatocellular carcinoma by promoting hepatocellular growth and survival (Tufan et al., 2002 [PubMed 12082552]).
Applications:	ELISA, IHC
Name of antibody:	URGCP
Immunogen:	Synthetic peptide of human URGCP
Full name:	upregulator of cell proliferation
Synonyms:	URG4
SwissProt:	Q8TCY9

ELISA Recommended dilution:	500-1000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

