

## 兔抗 NUDT12 多克隆抗体

中文名称：兔抗 NUDT12 多克隆抗体

英文名称：Anti-NUDT12 rabbit polyclonal antibody

相关类别：一抗

抗原：NUDT12

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Nucleotides are involved in numerous biochemical reactions and pathways within the cell as substrates, cofactors, and effectors. Nudix hydrolases, such as NUDT12, regulate the concentrations of individual nucleotides and of nucleotide ratios in response to changing circumstances.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	NUDT12
<b>Immunogen:</b>	Fusion protein of human NUDT12
<b>Full name:</b>	nudix (nucleoside diphosphate linked moiety X)-type motif 12
<b>SwissProt:</b>	Q9BQG2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	293T cell, K562 cell, Human pancreas tissue

WB Recommended dilution: 200-1000



