

NLGN1 抗原（重组蛋白）

中文名称： NLGN1 抗原（重组蛋白）

英文名称： NLGN1 Antigen (Recombinant Protein)

别 名： neurologin 1; NL1

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 44-300 amino acids of human NLGN1

技术规格

Full name:	neurologin 1
Synonyms:	NL1
Swissprot:	Q8N2Q7
Gene Accession:	BC032555
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.