

## 兔抗 NUDT18 多克隆抗体

中文名称：兔抗 NUDT18 多克隆抗体

英文名称：Anti-NUDT18 rabbit polyclonal antibody

别名：MTH3

相关类别：一抗

抗原：NUDT18

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the Nudix hydrolase family. Nudix hydrolases eliminate potentially toxic nucleotide metabolites from the cell and regulate the concentrations and availability of many different nucleotide substrates, cofactors, and signaling molecules. This protein contains a Nudix hydrolase domain and hydrolyzes oxidized forms of guanosine and deoxyguanosine diphosphates.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NUDT18
<b>Immunogen:</b>	Fusion protein of human NUDT18
<b>Full name:</b>	nudix (nucleoside diphosphate linked moiety X)-type motif 18

<b>Synonyms :</b>	MTH3
<b>SwissProt:</b>	Q6ZVK8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	30-150

