

NOC2L 抗原（重组蛋白）

中文名称： NOC2L 抗原（重组蛋白）

英文名称： NOC2L Antigen (Recombinant Protein)

别名： nucleolar complex associated 2 homolog; NIR; NET7; NET15

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 401-640 amino acids of human NOC2L

技术规格

Full name:	nucleolar complex associated 2 homolog
Synonyms:	NIR; NET7; NET15
Swissprot:	Q9Y3T9
Gene Accession:	BC003555
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Histone modification by histone acetyltransferases (HAT) and histone deacetylases (HDAC) can control major aspects of transcriptional regulation. NOC2L represents a novel HDAC-independent inhibitor of histone acetyltransferase (INHAT)