

兔抗 RPS6KA1 (Phospho-Thr348)多克隆抗体

- 中文名称：兔抗 RPS6KA1 (Phospho-Thr348)多克隆抗体
- 英文名称：Anti-RPS6KA1 (Phospho-Thr348) rabbit polyclonal antibody
- 别名：RSK; HU-1; RSK1; MAPKAPK1A
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：RPS6KA1 (Phospho-Thr348)
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signaling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the tr

	anscription factor CREB.
Applications:	WB, IHC, IF
Name of antibody:	RPS6KA1 (Phospho-Thr348)
Immunogen:	Synthetic peptide of human RPS6KA1 (Phospho-Thr348)
Full name:	ribosomal protein S6 kinase, 90kDa, polypeptide 1 (Phospho-Thr348)
Synonyms :	RSK; HU-1; RSK1; MAPKAPK1A
SwissProt:	Q15418
IHC positive control:	Human breast carcinoma
IHC Recommend dilution:	50-100
WB Predicted band size:	90 kDa
WB Positive control:	Hela cells untreated or treated with UV
WB Recommended dilution:	500-1000
IF Positive control:	Hela cells
IF Recommended dilution	100-200



