

兔抗 NUDT22 多克隆抗体

中文名称：兔抗 NUDT22 多克隆抗体

英文名称：Anti-NUDT22 rabbit polyclonal antibody

别名：nudix hydrolase 22

相关类别：一抗

储存：冷冻（-20℃）

抗原：NUDT22

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Hydrolyzes UDP-glucose to glucose 1-phosphate and UMP and UDP-galactose to galactose 1-phosphate and UMP . Preferred substrate is UDP-glucose.
Applications:	ELISA, WB
Name of antibody:	NUDT22
Immunogen:	Fusion protein of human NUDT22
Full name:	nudix hydrolase 22
SwissProt:	Q9BRQ3
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	33 kDa
WB Positive control:	231 cell, Jurkat cell, PC3 and 293T cell lysates
WB Recommended dilution:	1000-5000



