

兔抗 NUDT2 多克隆抗体

中文名称：兔抗 NUDT2 多克隆抗体

英文名称：Anti-NUDT2 rabbit polyclonal antibody

别名：APAH1

相关类别：一抗

抗原：NUDT2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the MutT family of nucleotide pyrophosphatases, a subset of the larger NUDIX hydrolase family. The gene product possesses a modification of the MutT sequence motif found in certain nucleotide pyrophosphatases. The enzyme asymmetrically hydrolyzes Ap4A to yield AMP and ATP and is responsible for maintaining the intracellular level of the dinucleotide Ap4A, the function of which has yet to be established. This gene may be a candidate tumor suppressor gene.
Applications:	ELISA, WB, IHC
Name of antibody:	NUDT2
Immunogen:	Full length fusion protein
Full name:	nudix (nucleoside diphosphate linked moiety X)-type moti

	f 2
Synonyms :	APAH1
SwissProt:	P50583
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	17 kDa
WB Positive control:	Jurkat cell, Mouse liver tissue and spleen tissue
WB Recommended dilution:	200-1000



