

## 兔抗 USP2 多克隆抗体

中文名称：兔抗 USP2 多克隆抗体

英文名称： Anti-USP2 rabbit polyclonal antibody

别名： USP9; UBP41

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： USP2

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a member of the family of de-ubiquitinating enzymes, which belongs to the peptidase C19 superfamily. The encoded protein is a ubiquitin-specific protease which is required for TNF-alpha (tumor necrosis factor alpha) -induced NF-kB (nuclear factor kappa B) signaling. This protein deubiquitinates polyubiquitinated target proteins such as fatty acid synthase, murine double minute 2 (MDM2), MDM4/MDMX and cyclin D1. MDM2 and MDM4 are negative regulators of the p53 tumor suppressor and cyclin D1 is required for cell cycle G1/S transition. Multip

	le alternatively spliced transcript variants encoding different isoforms have been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	USP2
<b>Immunogen:</b>	Fusion protein of human USP2
<b>Full name:</b>	ubiquitin specific peptidase 2
<b>Synonyms :</b>	USP9; UBP41
<b>SwissProt:</b>	O75604
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

