

兔抗 RPS6KA1 多克隆抗体

中文名称：兔抗 RPS6KA1 多克隆抗体

英文名称：Anti-RPS6KA1 rabbit polyclonal antibody

别名：RSK; HU-1; RSK1; MAPKAPK1A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RPS6KA1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Applications:	ELISA, WB, IHC
Name of antibody:	RPS6KA1

Immunogen:	Fusion protein of human RPS6KA1
Full name:	ribosomal protein S6 kinase, 90kDa, polypeptide 1
Synonyms :	RSK; HU-1; RSK1; MAPKAPK1A
SwissProt:	Q15418
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	83 kDa
WB Positive control:	K562 cells
WB Recommended dilution:	500-2000



