

## 兔抗 USP3 多克隆抗体

中文名称：兔抗 USP3 多克隆抗体

英文名称：Anti-USP3 rabbit polyclonal antibody

别名：UBP; SIH003

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：USP3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

USP3 (ubiquitin specific peptidase 3), also known as UBP or SIH003, is a 520 amino acid protein that contains one UBP-type zinc finger and belongs to the peptidase C19 family. Expressed ubiquitously with highest levels present in pancreas, USP3 catalyzes the conversion of a ubiquitin C-terminal thioester to a free ubiquitin and a thiol. The gene encoding USP3 maps to human chromosome 15, which houses over 700 genes and comprises nearly 3% of the human genome.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	USP3
<b>Immunogen:</b>	Fusion protein of human USP3
<b>Full name:</b>	ubiquitin specific peptidase 3
<b>Synonyms:</b>	UBP; SIH003
<b>SwissProt:</b>	Q9Y6I4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human gastric cancer
<b>IHC Recommend dilution:</b>	20-100



