

兔抗 RPS6KA2 多克隆抗体

中文名称：兔抗 RPS6KA2 多克隆抗体

英文名称： Anti-RPS6KA2 rabbit polyclonal antibody

别名： RSK; HU-2; RSK3; p90-RSK3; pp90RSK3; MAPKAPK1C; S6K-alpha; S6K-alpha2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPS6KA2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

IHC Recommend dilution:	50-200
IHC positive control:	Human liver cancer and Human colon cancer
ELISA Recommended dilution:	2000-5000
SwissProt:	Q15349
Synonyms :	RSK; HU-2; RSK3; p90-RSK3; pp90RSK3; MAPKAPK1C ; S6K-alpha; S6K-alpha2
Full name:	ribosomal protein S6 kinase, 90kDa, polypeptide 2
Immunogen:	Synthetic peptide of human RPS6KA2
Name of antibody:	RPS6KA2
Applications:	ELISA, IHC
Background:	This gene encodes a member of the RSK (ribosomal

S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signaling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

