

兔抗 ING2 多克隆抗体

中文名称：兔抗 ING2 多克隆抗体

英文名称：Anti-ING2 rabbit polyclonal antibody

别名：ING1L; p33ING2

抗原：ING2

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis.
Applications:	ELISA, IHC
Name of antibody:	ING2
Immunogen:	Fusion protein of human ING2
Full name:	inhibitor of growth family, member 2
Synonyms :	ING1L; p33ING2
SwissProt:	Q9H160

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human colon cancer
IHC Recommend dilution:	50-200

