

兔抗 RPS6KA4(Ab-568) 多克隆抗体

中文名称：兔抗 RPS6KA4(Ab-568) 多克隆抗体

英文名称： Anti-RPS6KA4(Ab-568) rabbit polyclonal antibody

别名： MSK2; RSK-B

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： RPS6KA4(Ab-568)

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates.

	ates, including CREB1 and c-fos. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.
Applications:	WB
Name of antibody:	RPS6KA4(Ab-568)
Immunogen:	Synthesized non-phosphopeptide derived from human MSK2 around the phosphorylation site of threonine 568 (M-Q-T(p)-P-C).
Full name:	ribosomal protein S6 kinase, 90kDa, polypeptide 4
Synonyms :	MSK2; RSK-B
SwissProt:	O75676
WB Predicted band size:	86 kDa
WB Positive control:	RAW264.7 cells treated with UV (5mins) lysate
WB Recommended dilution:	500-3000

