

## 兔抗 RPS6KA6 多克隆抗体

中文名称：兔抗 RPS6KA6 多克隆抗体

英文名称：Anti-RPS6KA6 rabbit polyclonal antibody

别名：RSK4; PP90RSK4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RPS6KA6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of ribosomal S6 kinase family, serine-threonine protein kinases which are regulated by growth factors. The encoded protein may be distinct from other members of this family, however, as studies suggest it is not growth factor dependent and may not participate in the same signaling pathways.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RPS6KA6
<b>Immunogen:</b>	Synthetic peptide of human RPS6KA6
<b>Full name:</b>	ribosomal protein S6 kinase, 90kDa, polypeptide 6

<b>Synonyms:</b>	RSK4; PP90RSK4
<b>SwissProt:</b>	Q9UK32
<b>IHC positive control:</b>	Human gastric cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

