

兔抗 USP5 多克隆抗体

中文名称: 兔抗 USP5 多克隆抗体

英文名称: Anti-USP5 rabbit polyclonal antibody

别名: ISOT

储存: 冷冻 (-20℃)

抗原: USP5

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Ubiquitin (see MIM 191339)-dependent proteolysis is a complex pathway of protein metabolism implicated in such diverse cellular functions as maintenance of chromatin structure, receptor function, and degradation of abnormal proteins. A late step of the process involves disassembly of the polyubiquitin chains on degraded proteins into ubiquitin monomers. USP5 disassembles branched polyubiquitin chains by a sequential exo mechanism, starting at the proximal end of the chain (Wilkinson et al., 1995 [PubMed 7578059]).
Applications:	ELISA, IHC

Name of antibody:	USP5
Immunogen:	Fusion protein of human USP5
Full name:	ubiquitin specific peptidase 5 (isopeptidase T)
Synonyms :	ISOT
SwissProt:	P45974
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100



