

NUP85 抗原（重组蛋白）

中文名称： NUP85 抗原（重组蛋白）

英文名称： NUP85 Antigen (Recombinant Protein)

别名： nucleoporin 85kDa; Nup75

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 250 amino acids of human NUP85

技术规格

Full name:	nucleoporin 85kDa
Synonyms:	Nup75
Swissprot:	Q9BW27
Gene Accession:	BC000697
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Bidirectional transport of macromolecules between the cytoplasm and nucleus occurs through nuclear pore complexes (NPCs) embedded in the nuclear envelope. NPCs are composed of subcomplexes, and NUP85 is part of one such subcomplex, Nup107-160.