

兔抗 USP6NL 多克隆抗体

中文名称：兔抗 USP6NL 多克隆抗体

英文名称：Anti-USP6NL rabbit polyclonal antibody

别名：RNTRE; TRE2NL; USP6NL-IT1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：USP6NL

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Acts as a GTPase-activating protein for RAB5A. Involved in receptor trafficking. In complex with EPS8 inhibits internalization of EGFR.
Applications:	ELISA, WB, IHC
Name of antibody:	USP6NL
Immunogen:	Fusion protein of human USP6NL
Full name:	USP6 N-terminal like
Synonyms:	RNTRE; TRE2NL; USP6NL-IT1
SwissProt:	Q92738
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100

WB Predicted band size:	94 kDa
WB Positive control:	293T cell lysate
WB Recommended dilution:	500-2000



